


RADIOLOGY DIGITAL IMAGING NETWORK PICTURE ARCHIVING AND COMMUNICATION SYSTEM (DINS-PACS)

The following are instructions on how to view x-rays VIA the DACH webserver:


1. Desktop PC must have version 5.5 or higher, Internet Explorer, 128 - bit, web browser
2. All Internet Explorer security settings located under **Tools** must be enabled. These changes only have the first time the user logs onto each PC. See below:



Click on Tools and select Internet Options
from the drop down menu




- 1. Select security**
- 2. Select Custom Level**




1. Scroll down and enable all security

2. Click OK when finished




- 1. Select Programs**
- 2. Check Internet Explorer as default web browser**
- 3. Reset web settings by clicking here**
- 4. Click OK when finished**
- 5. Proceed to Darnall Home Page at:
<http://143.83.0.101>**



Click on Radworks Webviewer




Click yes to security alert



- 1. If using window 2000 operating system enter your DIN-PACS Username, password, and **ppacs_gprmc** for domain**
- 2. If access is denied, call Mr Harris at 286-7650 or Ms Hall at 286-7167.**



**1. If using windows XP operating system enter
ppacs_gprmc\your username on the first line
your password on the second line**



1. To query by patient name enter patient's last name* followed by first name* (harris*jasper*) in name field

OR

2. To query by Patient ID enter the patient's full ID number (20/000-00-0000) in ID field

NOTE: Only query one field at a time

If study cannot be found by querying the

RadWorks WebViewer 5.1 - Microsoft Internet Explorer

File Edit View Favorites Tools Help

Back Forward Stop Refresh Home Search Favorites History Mail Size Print Edit

Address <https://dachradweb1.amedd.army.mil/ami/html/> Go

HARRIS, CARL, ANTONIO

11 of 11 records

DACH Webserver

Series

Series

☐ Key Images only

View Selection

Quality

No compression

Loc.	Patient	Access.Nun
+	HARRIS, ANGELA, M	0110/010
+	HARRIS, CARL, ANTONIO	0110/010
+	HARRIS, CHARMAYNE, ALANA	0110/010
+	HARRIS, CHARMAYNE, ALANA	0110/010
+	HARRIS, CHARMAYNE, ALANA	0110/010
+	HARRIS, GERALD, N	0110/010
+	HARRIS, JASON, LAMONT	0110/010
+	HARRIS, JIMMY, EDWARD	0110/010

Modality

☒ CR

☒ NM

☒ US


☒ unknown

All None


Report DACH Home

Done Internet

1. To view selected study double click on the the patient's name



1. To magnify, flip, and/or rotate image left click on this tool tab
2. To view CT or MRI studies in cine left click on this tool tab
3. To measure distance left click on this tool tab
4. To change how images are displayed left click on this tool tab



Menu choices to zoom, magnify, flip/rotate, measure, etc can also be accessed by a right of the mouse any where on the image.